

## How to install LANDesk on Windows OS

1. Connect to the MSU network via a wired connection or MSU-Secure wireless.
2. Install the LANDesk agent by running this package:  
`\\ldms\agents\MSUWinbase_with_status.exe`
3. Reboot (recommended) after installation is complete.

## How to install LANDesk on Mac OS X

Before you begin you must be connected to MSU network via a wired connection or MSU-Secure wireless.

1. Open **Finder** and click on **Go** menu.
2. Select **Connect to Server...**
3. Enter `smb://ldms.msu.montana.edu` in the Server Address field and click the **Connect** button.
4. *If prompted for credentials, enter your NetID and password*
5. Select **agents** in the “Select the volumes you want to mount...” window and then click the **OK** button.
6. Drag or copy the `LDMSClient.mpkg` to your desktop.
7. *Wait for download to complete before going to step 6.*
8. Click the `LDMSClient.mpkg` to open and then click the **Continue** button.
9. Click the **Continue** button to select the destination drive, generally Macintosh HD.
10. Click the **Install** button. Enter your computer’s credentials when prompted.
11. Click the **Close** button (after the installation has been successful).
12. *The `LDMSClient.mpkg` can now be moved to the trash.*
13. Click the **LANDesk icon** in the upper toolbar on the screen and select **Inventory Scan**.

Need Help? Please contact the IT Service Desk (994-1777 | [helpdesk@montana.edu](mailto:helpdesk@montana.edu)).